

Medical Outlet
MGS-CL 301




DESCRIPTION

Gas sockets and jacks are produced in BS, DIN, NF / AF norms. Jacks can be used with suitable standard sockets. Sockets and Jacks are produced according to EN ISO 9170-1.

STRUCTURAL FEATURES

Dimensions	~5,5x2 cm
Parts	Base and Head
Weight	~0,3Kg
Norm	DIN, BS, AF
Oxygen(O ₂)	5bar DIN - 3,5bar AF - 4barBS
Vacuum(VAC)	-0,9bar
Air(H ₄ ,H ₇)	5bar DIN - 3,5bar AF - 4barBS
NitrousOxide(N ₂ O)	40 ±3,1/min
Protection Class	Class I
Medical Device Class	Class 2b Rule11
Colour	
Hose Inlets	4 mm, 6 mm
Dimensions of box	70mm x 90mm x 360mm



Medical Gas Outlets



PRODUCT CODES

SUBMODELS REFERENCE NO	MODEL NAME	SUBMODELS REFERENCE NO	MODEL NAME
250-01-01	MEDICAL GAS OUTLET O2 DIN	250-08- 01	MEDICAL GAS OUTLET N2 DIN
250-01-02	MEDICAL GAS OUTLET O2 BS	250-08- 02	MEDICAL GAS OUTLET N2 BS
250-01-03	MEDICAL GAS OUTLET O2 AF	250-08- 03	MEDICAL GAS OUTLET N2 AF
250-01-04	MEDICAL GAS OUTLET O2 OHMEDA	250-08- 04	MEDICAL GAS OUTLET N2 OHMEDA
250-01-05	MEDICAL GAS OUTLET W/ UNION OUTLET O2 DIN	250-08- 05	MEDICAL GAS OUTLET W/ UNION OUTLET N2 DIN
250-01-06	MEDICAL GAS OUTLET W/ UNION OUTLET O2 BS	250-08- 06	MEDICAL GAS OUTLET W/ UNION OUTLET N2 BS
250-01-07	MEDICAL GAS OUTLET W/ UNION OUTLET O2 AF	250-08- 07	MEDICAL GAS OUTLET W/ UNION OUTLET N2 AF
250-01-08	MEDICAL GAS OUTLET W/ UNION OUTLET O2 OHMEDA	250-08- 08	MEDICAL GAS OUTLET W/ UNION OUTLET N2 OHMEDA
250-01-09	MEDICAL GAS OUTLET PNEUMATIC O2 DIN	250-08- 09	MEDICAL GAS OUTLET PNEUMATIC N2 DIN
250-01-10	MEDICAL GAS OUTLET PNEUMATIC O2 BS	250-08- 10	MEDICAL GAS OUTLET PNEUMATIC N2 BS
250-01-11	MEDICAL GAS OUTLET PNEUMATIC O2 AF	250-08- 11	MEDICAL GAS OUTLET PNEUMATIC N2 AF
250-01-12	MEDICAL GAS OUTLET PNEUMATIC O2 OHMEDA	250-08- 12	MEDICAL GAS OUTLET PNEUMATIC N2 OHMEDA
250-02-01	MEDICAL GAS OUTLET N2O DIN	250-09- 01	MEDICAL GAS OUTLET DUPLICATOR O2 DIN
250-02-02	MEDICAL GAS OUTLET N2O BS	250-09- 02	MEDICAL GAS OUTLET DUPLICATOR O2 BS
250-02-03	MEDICAL GAS OUTLET N2O AF	250-09- 03	MEDICAL GAS OUTLET DUPLICATOR O2 AF
250-02-04	MEDICAL GAS OUTLET N2O OHMEDA	250-09- 04	MEDICAL GAS OUTLET DUPLICATOR O2 OHMEDA
250-02-05	MEDICAL GAS OUTLET W/ UNION OUTLET N2O DIN	250-10- 01	MEDICAL GAS OUTLET DUPLICATOR N2O DIN
250-02-06	MEDICAL GAS OUTLET W/ UNION OUTLET N2O BS	250-10- 02	MEDICAL GAS OUTLET DUPLICATOR N2O BS
250-02-07	MEDICAL GAS OUTLET W/ UNION OUTLET N2O AF	250-10- 03	MEDICAL GAS OUTLET DUPLICATOR N2O AF
250-02-08	MEDICAL GAS OUTLET W/ UNION OUTLET N2O OHMEDA	250-10- 04	MEDICAL GAS OUTLET DUPLICATOR N2O OHMEDA
250-02-09	MEDICAL GAS OUTLET PNEUMATIC N2O DIN	250-11- 01	MEDICAL GAS OUTLET DUPLICATOR MA4 DIN
250-02-10	MEDICAL GAS OUTLET PNEUMATIC N2O BS	250-11- 02	MEDICAL GAS OUTLET DUPLICATOR MA4 BS
250-02-11	MEDICAL GAS OUTLET PNEUMATIC N2O AF	250-11- 03	MEDICAL GAS OUTLET DUPLICATOR MA4 AF
250-02-12	MEDICAL GAS OUTLET PNEUMATIC N2O OHMEDA	250-11- 04	MEDICAL GAS OUTLET DUPLICATOR MA4 OHMEDA
250-03-01	MEDICAL GAS OUTLET MA4 DIN	250-12- 01	MEDICAL GAS OUTLET DUPLICATOR SA7 DIN
250-03-02	MEDICAL GAS OUTLET MA4 BS	250-12- 02	MEDICAL GAS OUTLET DUPLICATOR SA7 BS
250-03-03	MEDICAL GAS OUTLET MA4 AF	250-12- 03	MEDICAL GAS OUTLET DUPLICATOR SA7 AF
250-03-04	MEDICAL GAS OUTLET MA4 OHMEDA	250-12- 04	MEDICAL GAS OUTLET DUPLICATOR SA7 OHMEDA



PRODUCT CODES



700-01-01	HOSE OUTLET JACK (PROBE) O2 DIN
700-01-02	HOSE OUTLET JACK (PROBE) O2 BS
700-01-03	HOSE OUTLET JACK (PROBE) O2 AF
700-01-04	HOSE OUTLET JACK (PROBE) O2 OHMEDA
700-02-01	HOSE OUTLET JACK (PROBE) N2O DIN
700-02-02	HOSE OUTLET JACK (PROBE) N2O BS
700-02-03	HOSE OUTLET JACK (PROBE) N2O AF
700-02-04	HOSE OUTLET JACK (PROBE) N2O OHMEDA
700-03-01	HOSE OUTLET JACK (PROBE) MA4 DIN
700-03-02	HOSE OUTLET JACK (PROBE) MA4 BS
700-03-03	HOSE OUTLET JACK (PROBE) MA4 OHMEDA
700-03-04	HOSE OUTLET JACK (PROBE) MA4 AF
700-04-01	HOSE OUTLET JACK (PROBE) SA7 DIN
700-04-02	HOSE OUTLET JACK (PROBE) SA7 BS
700-04-03	HOSE OUTLET JACK (PROBE) SA7 AF
700-04-04	HOSE OUTLET JACK (PROBE) SA7 OHMEDA
700-07-01	HOSE OUTLET JACK (PROBE) CO2 DIN
700-07-02	HOSE OUTLET JACK (PROBE) CO2 BS
700-07-03	HOSE OUTLET JACK (PROBE) CO2 AF
700-07-04	HOSE OUTLET JACK (PROBE) CO2 OHMEDA
700-08-01	HOSE OUTLET JACK (PROBE) N2 DIN
700-08-02	HOSE OUTLET JACK (PROBE) N2 BS
700-08-03	HOSE OUTLET JACK (PROBE) N2 AF
700-08-04	HOSE OUTLET JACK (PROBE) N2 OHMEDA
700-01-01	HOSE OUTLET JACK (PROBE) O2 DIN
700-01-02	HOSE OUTLET JACK (PROBE) O2 BS
700-01-03	HOSE OUTLET JACK (PROBE) O2 AF
700-01-04	HOSE OUTLET JACK (PROBE) O2 OHMEDA
700-02-01	HOSE OUTLET JACK (PROBE) N2O DIN
700-02-02	HOSE OUTLET JACK (PROBE) N2O BS
700-02-03	HOSE OUTLET JACK (PROBE) N2O AF
700-02-04	HOSE OUTLET JACK (PROBE) N2O OHMEDA
700-03-01	HOSE OUTLET JACK (PROBE) MA4 DIN
700-03-02	HOSE OUTLET JACK (PROBE) MA4 BS
700-03-03	HOSE OUTLET JACK (PROBE) MA4 OHMEDA
700-03-04	HOSE OUTLET JACK (PROBE) MA4 AF
700-04-01	HOSE OUTLET JACK (PROBE) SA7 DIN
700-04-02	HOSE OUTLET JACK (PROBE) SA7 BS
700-04-03	HOSE OUTLET JACK (PROBE) SA7 AF
700-04-04	HOSE OUTLET JACK (PROBE) SA7 OHMEDA
700-07-01	HOSE OUTLET JACK (PROBE) CO2 DIN
700-07-02	HOSE OUTLET JACK (PROBE) CO2 BS
700-07-03	HOSE OUTLET JACK (PROBE) CO2 AF
700-07-04	HOSE OUTLET JACK (PROBE) CO2 OHMEDA
700-08-01	HOSE OUTLET JACK (PROBE) N2 DIN
700-08-02	HOSE OUTLET JACK (PROBE) N2 BS
700-08-03	HOSE OUTLET JACK (PROBE) N2 AF
700-08-04	HOSE OUTLET JACK (PROBE) N2 OHMEDA



+44 7401 868130



info@cliniqonhc.uk
www.cliniqonhc.uk



Hamilton House,
4a The Avenue,
London, England, E4 9LD