CLINIQON HC



LEADING INNOVATIONS



Mutipurpose Chairs
GRACE MPC 250

GRACE MPC 250

PRODUCT FEATURES

- Back, foot and lowering-uppering height adjustment are done by 4 pieces of DC motors and the desired position is provided.
- Mattress is made by viscoelastic foam and is designed to take the shape of patients body
 when the seating for long time
- Head, back, seat, foot and arm parts of the upholstery can be taken down and put on easily to clean
- The arms consist of 2 parts, the elbow parts which are fixed while the front parts are able to move to the left – right, up – down within 3 steps according to the patient position during the dialyzing process.
- The arms, related to the back angle, can move without disturbing their parallel to the ground.
- There is a foot rest in foot part of the chair as taking down-putting on the type to prevent the cramps of the patients during the dialysis procedure.
- The chair has diagonal brake system, designed to be locked when needed.
- The height of the chair head can be adjusted to the position of the patient, and it can be removable

• The steel construction is protected due to corrosion damage with the galvanize covering and electrostatic powder paint and is covered by antibacterial, easy cleaning feature ABS plastics.

• Back Segment Length: 77 cm

• Back Segment Angle: 0-75°

Min-Max Height: 59-79 cm

• Foot Segment Angle: 0-35°

• There are 3 parts as seat, back and foot.

• Seat Length: 50 cm

• Foot Segment Length: 64,5 cm

• Foot Support: 27,5 cm

• Width: 99 cm

Length: 205 cm



CLINIQON HC



LEADING INNOVATIONS



Mutipurpose Chairs
GRACE MPC 250 M

GRACE MPC 250 M

PRODUCT FEATURES

• Overall Width: 76 cm

• Height: 60 - 80 cm

• Overall Length: 115 cm

• Seating Surface Length: 185 cm

• Seating Surface Width: 60 cm

• 4 motors with height, backrest, foot and Trendelenburg adjustable

• Steel frame

• Electrostatic painted frame

• Adjustable with hand control

• Antibacterial, wipeable leather upholstery

Removable pillow for face opening







